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EMS

Bibliografie

1. Axelsson P, Lindhe J: The effect of a preventive programme on dental plaque, gingivitis and caries on schoolchildren. *J Clin Periodontol* 1, 126-138 (1974).
2. Axelsson P, Lindhe J: Effect on controlled oral hygiene procedures on caries and periodontal disease in adults. *J Clin Periodontol* 5, 133-151 (1978).
3. Axelsson P, Lindhe J: Effect of controlled oral hygiene procedures on caries and periodontal disease in adults. Results after 6 years. *J Clin Periodontol* 8, 239-248 (1981).
4. Bastendorf K-D, Bartsch A: Langzeiterfolge der Prophylaxe nach 30 Jahren Recall. *Prophylaxe impuls* 16, 62-69 (2012).
5. Strafela-Bastendorf N, Bastendorf KD: PZR-neu gedacht! *zm* 106 (11A), 26-30 (2016).
6. Strafela-Bastendorf N, Bastendorf KD: Professionelle Biofilamentfernung-Tipps für den Praxisalltag. *Quintessenz* 67, 1069-1075 (2016).
7. Wennström JL, Tomasi C, Bertelle A, Dellasega E: Full mouth ultrasonic debridement versus quadrant scaling and root planing as an initial approach in the treatment of chronic periodontitis *J Clin Periodontol* 32, 851-859 (2005). doi:10.1111/j.1600-051X.2005.00776.x.
8. M. Aslund, Suvan J, Moles DR, Dàiuo F, Tonetti MS: Effects of Two Different Methods of Non-Surgical Periodontal Therapy on Patient Perception of Pain and Quality of Life: A Randomized Controlled Clinical Trial. *J Periodontol* (2008). doi: 10.1902/jop.2008.070394
9. Wennström JL, Dahlen G, Ramberg P: Subgingival debridement of periodontal pockets by air polishing in comparison with ultrasonic instrumentation during maintenance therapy. *J Clin Periodontol* 38, 820-827 (2011). doi:10.1111/j.1600-051X.
10. Bühler J, Amato M, Weiger R, Walter C.: A systematic review on the effects of airpolishing devices on oral tissues. *Int J Dent Hygiene* (2015). doi:10.1111/idh.12120.
11. Sultan DA, Hill RG, Gillam DG: Air-Polishing in Subgingival Root Debridement: A Critical Literature Review. *Journal of Dentistry and Oral Biology* 2 (10), Article 1065 (2017).
12. Ethan N, Roy B, Spahr A, Tihana DR: The efficacy of air polishing devices in supportive periodontal therapy: A systematic review and meta-analysis. *Quintessenz International* (2018). doi:10.3290/j.qi.a.40341.
13. Moëne R, Décaillet F, Andersen E, Mombelli A: Subgingival plaque removal using a new air-polishing device. *Journal of Periodontology* 81, 79-88 (2010).
14. Kim MJ et al: Efficiency of professional tooth brushing before ultrasonic scaling (2015). doi: 10.1111/idh.12127.
15. Baehni, P.: Anwendung von Mundspülungen im Dentalbereich, *Prophylaxedialog* 2/1, 17-22 (2008/2009).
16. Sawhney A, Venugopol S, Babu G, Garg A, Mathew M, Yadav M, Gupta B, Tripathi S: Aerosols how dangerous they are in clinical practice. doi: 10.7860/JCDR/2015/12038.5835.
17. Bastendorf-Strafela N, Bastendorf KD, Mann P: Pilotstudie: Verbessert das Anfärben der Plaque die Ergebnisse der PZR? *PNC* (10) 2, 91-93 (2016).
18. Volgenant CMC, Fernandez y Mostajo M, Rosema NAM et al. Comparison of red autofluorescing plaque and disclosed plaque – a cross-sectional study. *Clin Oral Invest* 20, 2551-2558 (2016).
19. Mensi M, Agosti R, Cappa V, Calza S: The efficacy of disclosing plaque agent as a guide to the supra-gingival biofilm removal. Poster EURO PERIO in London (2014).
20. Schiffner U: Die Einwirkungen eines abrasiven Pulverstrahles auf artifiziell demineralisierten Schmelz. *Dtsch Zahnärztl Z* 47, 778-781 (1999).
21. Petersilka GJ et al (2003b): Subgingival plaque removal in buccal and lingual sites using a novel low abrasive air-polishing powder. *Journal of Clinical Periodontology* 30, 328-333 (2003).
22. Petersilka GJ et al (2003a): Subgingival plaque removal at interdental sites using a low abrasive air polishing powder. *Journal of Periodontology* 74, 307-311 (2003).
23. Flemmig TF et al: Subgingival debridement efficacy of glycine powder air-polishing. *Journal of Periodontology* 78, 1002-1010 (2007).
24. Flemmig TF et al: Randomized controlled trial assessing efficacy and safety of glycine powder air-polishing in moderate to deep periodontal pockets. *Journal of Periodontology* 83, 444-452 (2012).
25. Hägi T, Hofmänner P, Eick S, Donnet M, Salvi G, Sculean A, Ramseier C: The effects of erythritol air-polishing powder on microbiological and clinical outcomes during supportive periodontal therapy. Six months results of a randomized controlled clinical trial. *Quintessenz Int* 46, 31-41 (2015).
26. Sculean A, Bastendorf K.-D. et al: A paradigm shift in mechanical biofilm management? Subgingival air polishing: a new way to improve mechanical biofilm management in the dental practice. *Quintessenz International* 44 (7), (2013).
27. Graetz C et al: Removal of simulated biofilm: an evolution of the effect on root surfaces roughness after scaling. *Clin Oral Invest* (2016). doi: 10.1007/s00784-016-1861-9.
28. Busslinger A, Lampe K, Beuchat M, Lehmann B: A comparative in vitro study of the magnetostruktive und piezoelektrische ultrasonische Scaling Instrumente. *J Clin Periodontol* 28, 642-649 (2001).
29. Barendregt DS et al: Penetration depths with an ultrasonic mini insert compared with a conventional curette in patients with periodontitis and in periodontal maintenance. *J Clin Periodontol* 35, 31-36 (2008).
30. Flemmig TF, Petersilka GJ, Mehl A, Hickel R, Klaiber B: The effect of working parameters on root substance removal using a piezoelektrische ultrasonische Scaler in vitro. *J Periodontol* 25, 158-163 (1998a).
31. Flemmig TF, Petersilka GJ, Mehl A, Rudiger S, Hickel R, Klaiber B: Working parameters of a sonic scaler influencing root substance removal in vitro. *Clin Oral Invest* 1, 55-60 (1997).
32. Flemmig TF, Petersilka GJ, Mehl A, Hickel R, Klaiber B: Working parameters of a magnetostruktive ultrasonische Scaler influencing root surface removal in vitro. *J Periodontol* 69, 547-553 (1998b).

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© EMS